

Strategies for the Management of Microbial Contamination in Laboratory Incubators

Introduction:

Microbial contamination poses a significant challenge in cell culture laboratories, compromising research integrity and leading to inaccurate experimental results. For example, estimates for the incidence of *Mycoplasma* sp. contamination vary from 10% to 36% of cell lines used in laboratory procedures. Laboratory incubators, being a critical environment for cell growth, are prone to mycoplasma contamination. This paper explores effective strategies to eliminate mycoplasma contamination from laboratory incubators, ensuring reliable and consistent cell culture experiments.

I. Cleaning and Decontamination Procedures:

To minimize microbial contamination, thorough cleaning and decontamination protocols are essential. The following steps can be employed:

1. **Regular Cleaning:** Implement a stringent cleaning regimen for laboratory incubators, including the removal of all cell culture vessels, and cleaning of shelves, trays and support hardware. Clean surfaces using detergents (laboratory detergents such as Alconox or equivalent) to eliminate gross dirt and grime.
2. **Disinfection:** After cleaning, disinfect the incubator using (eg) heat treatment, hydrogen peroxide vapor (nebulized or vaporized H₂O₂), paracetic acid (PAA), or ultraviolet (UV) light. These methods have varying degrees of efficacy so a lab should choose based their own needs.
3. **Periodic Maintenance:** Establish a maintenance schedule for incubators to ensure proper functioning and prevent mycoplasma contamination. This includes regular replacement of filters, inspection of seals, and calibration of temperature and humidity controls.

II. Isolation and Quarantine:

Preventing microbial contamination requires isolating potentially contaminated cell lines and practicing strict quarantine measures:

1. **Isolation:** Identify and segregate contaminated cell lines from healthy ones. Store contaminated cultures separately, preferably in sealed containers, to prevent cross-contamination.
2. **Quarantine:** Newly obtained or suspected cell lines should be quarantined for mycoplasma testing before integration into the laboratory. Maintain a dedicated space



for quarantine and perform regular testing to ensure early detection and prompt elimination of contamination sources.

3. In some cases it may be possible to recover *Mycoplasma*-contaminated cultures using fluoroquinolone antibiotics. If there is no other source for uncontaminated cells, the time and expense may be justified.
4. Destroy remaining contaminated cultures by autoclaving.

III. Regular Testing:

In some cases, such as *Mycoplasma* contamination, rigorous testing protocols are indispensable in identifying mycoplasma contamination and taking appropriate measures:

1. Testing Methods: Employ validated detection methods, such as contact or settle plates or swabbing, polymerase chain reaction (PCR) or DNA staining assays. These tests allow for accurate detection and identification of contamination.
2. Frequency: If a contamination is detected, immediate attempts at isolation and identification should be made.
3. External Testing: Periodically utilize external services or contract laboratories for recurring or problematic contaminations.

Conclusion:

Microbial contamination represents a persistent challenge in laboratory incubators and can compromise the validity of cell culture experiments. By implementing robust cleaning and decontamination procedures, ensuring isolation and quarantine practices, and implementing regular microbial testing, laboratories can significantly reduce the risk of contamination. Maintaining a vigilant approach towards microbial elimination is crucial for generating reliable and reproducible experimental results, thus promoting the advancement of scientific research.

References

- Drexler H, Uphof CC (2002). *In Vitro Cell. Biol. - Animal* 38:79-85
- Olarerin-George AO, Hogenesch JB (2015). *Nucleic Acids Research* 43(5): 2535-2542
- Uphof CC, Denkman Sabine-A, Drexler HG (2012) *J. Biomed. Biotechnol.* 2012: 267678
- Chat-GPT 4, OpenAI

